

# Next Generation TruClip Holder

A proficient solution to your clinical workflow in pressure monitoring.



The next generation Edwards TruClip holder provides a streamlined and efficient sensor organization that can be used across the continuum of care – from OR into the ICU.


## Simple, convenient and easy to use

The Edwards Lifesciences TruClip holder quickly attaches to a wide variety of IV poles and can hold up to seven Edwards pressure monitoring, arterial pressure cardiac output (APCO) and blood sampling devices.

The TruClip holder's design facilitates quick setup by replacing traditional C-clamp holders that require two hands to operate. A simple squeeze-and-release clamping action allows single-handed operation for efficient mounting and quick adjustments in time-sensitive situations. The reusable TruClip holder is easy to use, has a reliable grip, and is adaptable across a range of departments.

### Now offering:

- Color-coded inserts
- Improved rubber grips
- Design that fits a wider variety of IV poles
- Four mounting spaces on the front



### Efficient sensor organization.

The TruClip holder includes color-coded inserts for easy identification of disposable pressure transducer (DPT) lines.



### Holds up to seven Edwards devices.

The TruClip holder is configured for the efficient organization of pressure monitoring and blood sampling devices. Four mounting spaces on the front help provide a clear, uncluttered environment in the OR and at the bedside.

### Reliable

Color-coded inserts for easy identification of DPT lines

### Durable

Improved rubber grips, reinforced structure, robust design

### Compatible

Accommodates a wide variety of IV poles

### Ease of use

Intuitive design offers quick and easy attachment and adjustment using just one hand



Edwards

## Variety of use environments

The TruClip holder can be used to help reduce equipment clutter in the OR, ICU, ED and Cardiac Catheter labs.



TruClip holder shown with Edwards VAMP closed blood sampling system.

## Maximum efficiency with your clinical workflow

By optimizing time to allow quick and simple organization of bedside systems, the TruClip holder may improve the efficiency of invasive monitoring throughout the hospital.

The versatile design allows the TruClip holder to work with a variety of devices, including the Edwards TruWave disposable pressure transducers, VAMP Adult and VAMP Plus closed blood sampling systems, Z-site sampling ports, FloTrac sensors and VolumeView sensors.

Model number	Dimensions	Weight limit	Compatible IV pole diameters
TCLIP05	8" L x 2.75" W x 2.5" H	2.5 lbs	0.5" – 1.5"

**For over 40 years, Edwards Lifesciences has been helping you make proactive clinical decisions to advance the care of surgical and critical care patients.**

Through ongoing collaboration with you, ongoing education and our never-ending quest for advancement, Edwards develops solutions that provide the clarity to make proactive clinical decisions.

**Know more.** [Edwards.com/TruClip](http://Edwards.com/TruClip)

**For professional use. CAUTION: Federal (United States) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.**

Edwards Lifesciences devices placed on the European market meeting the essential requirements referred to in Article 3 of the Medical Device Directive 93/42/EEC bear the CE marking of conformity.

Edwards, Edwards Lifesciences, the stylized E logo, FloTrac, TruClip, TruWave, VAMP, VAMP Plus, and VolumeView are trademarks of Edwards Lifesciences Corporation or its affiliates. All other trademarks are the property of their respective owners.

© 2016 Edwards Lifesciences Corporation. All rights reserved. PP--US-1766 v1.0

Edwards Lifesciences • One Edwards Way, Irvine CA 92614 USA • [edwards.com](http://edwards.com)

